

Stock Code: 300298

For Clinic

Passed the Domestic First-class
Hospital Clinical Test

Sinocare

• caring for love

Safe AQ Max I

Blood Glucose Monitor



5s Test Time



Easy to Use



Ketone Warning



500 Results Memory



Wide HCT Range:10%~70%



Match EN ISO 15197:2015

Changsha Sinocare Inc.



NMPA
(CFDA)

CE 0123

Safe AQ Max I

Blood Glucose Monitor



SPECIFICATION

FAD-GDH System

Blood Sample: Capillary blood, Venous blood,
Arterial blood, Neonatal blood

Calibration Sample: Venous blood

Test Time: 5s

Blood Volume: 1 μ L

HCT: 10%~70%

Test Range: 20~600 mg/dL, 1.1~33.3 mmol/L

Weight: 53.8g without battery

Dimensions: 85mm(L) x53mm(W)x23.6mm(H)

Memory: 500 blood glucose test results, 100 control solution test results



1 Insert Safe AQ Max I
test strip



2 Apply the blood sample



3 Result will be displayed in
5 seconds

Changsha Sinocare Inc.

No.265, Guyuan Road, Hi-tech Zone, Changsha, Hunan Province, 410205, People's Republic of China
TEL:+86-731-89935581/89935582 Email:info@sinocare.com Website:www.sinocare.com

Sinocare

• Caring for love