

Stock Code: 300298

**Sinocare**  
caring for love

# Safe AQ Max II

Blood Glucose Monitoring System



 **Bluetooth**



## Smart

Color Indicator, Beep Function  
Helpful Meal Markers, Ketone Warning



## Accurate

Proven Accuracy Match EN ISO 15197:2015  
16 Sliver & Carton Electrodes System  
New Fully-automatic Laser Correction Code Technology



## Fast

5s Test Time



## Easy

Convenient Auto Strip & Lancet Ejector  
Large Display

Changsha Sinocare Inc.



NMPA  
(CFDA)

CE 0123



# Safe AQ Max II

## Blood Glucose Monitoring System



### SPECIFICATION

FAD – GDH Enzyme

Blood Sample: Capillary blood, Venous blood,  
Arterial blood, Neonatal blood

Measurement units: mg/dL or mmol/L

Measurement range: 20~600 mg/dL (1.1~33.3 mmol/L)

Wide 10%-70% Hematocrit(HCT)

Tiny 1 µl blood sample

Fast 5 seconds countdown

No Coding

7-, 14-, 30-, and 90-day averaging

Ketone Warning

Lightning Indicator

Event Tagging: after meal/ before meal/ random

Stores 500 test results and 90 control results

Test Temperature: 5°C~45°C

Test related humidity: ≤90%

On/Off Sound

2\* AAA Alkaline Batteries



1 Insert Safe AQ Max II test strip



2 Collect blood sample



3 Apply the blood sample



4 Test results displayed in 5 seconds



5 Discard the used test strip & lancet



Changsha Sinocare Inc.

No.265, Guyuan Road, Hi-tech Zone, Changsha, Hunan Province, 410205, People's Republic of China  
TEL:+86-731-89935581/89935582 Email:info@sinocare.com Website:www.sinocareintl.com

**Sinocare**

• Caring for love