

Stock Code: 300298

Sinocare
caring for love

Safe AQ UG

Blood Glucose & Uric Acid Monitoring System



- Dual function
- Tiny sample volume
- Smooth strip ejector
- Anti-interference
- Clinically proven accuracy

Sinocare Inc.



CE 0123

Safe AQ UG

Blood Glucose & Uric Acid Monitoring System

SPECIFICATION

Sample type : Capillary or venous whole blood

Sample volume : blood glucose 0.6 μ L, uric acid 3 μ L


Calibration : Plasma equivalent

Measuring range : Blood glucose 20-600mg/dL
or 1.1-33.3mmol/L, uric acid
3.0-20.0mg/dL or 181-1188 μ mol/L

Test time : Blood glucose 5 $\pm$ 1seconds,
uric acid 25 $\pm$ 1seconds

Memory capacity : Blood glucose 500 test results and 100 control solution
with date and time, uric acid 100 test results and 100
control solution with date and time

Battery life : 1000 tests




① Insert uric acid test strip



② Apply the blood sample



③ Result will be displayed in 25 seconds



① Insert blood glucose test strip



② Apply the blood sample



③ Result will be displayed in 5 seconds

Sinocare Inc.

No.265, Guyuan Road, Hi-tech Zone, Changsha, Hunan Province, 410205, People's Republic of China
TEL: +86-731-89935581/89935582 Email: info@sinocare.com Website: www.sinocareintl.com

Sinocare
• Caring for love